

## ABDOS Elegant Pipette: light, safe & reliable

### Fully autoclavable micropipette for everyday laboratory use

The Elegant pipette from our partner ABDOS is a fully **autoclavable**, ergonomically designed micropipette that offers a high level of **accuracy** and **precision**. The pipette is ultra-light and **requires** only a **low pipetting force**. The soft volume setting and smooth plunger movement allow effortless and safe pipetting. Each pipette is supplied with an individual CE test certificate and is calibrated according to ISO 8655 in conjunction with ISO 13485.



The **Elegant Pipette** is available in two versions:


- Volume **2 – 20 µl**, adjustable in 0,2 µL-steps
- Volume **20 – 200 µl**, adjustable in 2,0 µL-steps

Both pipettes can be calibrated using a tool supplied and can be easily distinguished by different colored caps.

Find out more: <https://www.biomed.de/en/products/laboratory-supplies/>

# Accessories: pipette tips & pipette holders


## Pipette tips in a useful & recyclable holder



The practical and reusable **pipette tip holder** with the **CE & IVD certified** pipette tips are already inserted and make everyday laboratory work easier. The holders are neatly stored in a box that can be **autoclaved** up to ten times. The pipette tips can be securely attached to the Elegant Pipette or other commercially available micropipettes. They seal tightly and are absolutely leak-proof. In addition, the pipette tips are equipped with a special filter that **prevents contamination**.

They are made of **high purity polypropylene (PP)** and are **hydrophobic** according to US FDA 21CFR. They are **free** of **DNase** and **RNase**, **endotoxins**, **ATP & pyrogens**, and **human DNA**, and are also free of **heavy metals**, **natural rubber**, and **inhibitory plasticizers**, as well as **RoHS-free** according to EU Directive 2011/65/EU.

The pipette tips and clear box comply with ISO 9001:2015 in conjunction with ISO 13485:2016, and are available in **2 - 20 µl** and **2 - 200 µl** sizes.



## For Organization & Cleanliness in the Lab: Pipette Stand & Carousel

The practical pipette rack or the space-saving pipette carousel provide a tidy workplace.



### Pipette rack

The practical pipette stand made of acrylic glass offers space for **up to five pipettes**. It is easy to clean and has a safe and firm stand.

### Pipette carousel

Allow your pipettes to have some fun after work with the rotating pipette carousel. **Up to six pipettes** can be stowed in the space-saving pipette carousel. The carousel has a firm stand and can be rotated smoothly and evenly. It is also **autoclavable**.



Visit our website: <https://www.biomed.de/en/products/laboratory-supplies/>